




PRODUCT CATALOG 2018

Enabling
the Future
through
Nanotechnology



www.nanointegris.com
sales@nanointegris.com
T: 1.866.650.0482





2018 Product Catalog

Brief History

In October 2006, Professor Mark Hersam's research group at Northwestern University published a ground breaking paper in Nature Nanotechnology describing a process to sort CNTs by electronic structure. Flooded with sample requests from around the world, NanoIntegris was founded in January 2007 and established in Skokie, Illinois. Raymor, a high-value added materials supplier, acquired NanoIntegris in 2012 to expand its client base and its expertise in nanotube processing.

Over the past 10 years, Raymor Industries has developed its plasma processing capability which led to the marketing of two outstanding products: plasma-grown single-wall carbon nanotubes (SWCNT) and plasma-grown crumpled graphene nanoplatelets.

Raymor and NanoIntegris offer a wide variety of premium nanomaterials to companies and academic institutions developing next-generation electronics, energy, and biomedical technologies. We pride ourselves on our thorough, accurate, and honest material characterization.

Our nanotube powders, inks for printed electronics, and dispersions are among the purest in the industry. What is more, our strict quality control standards and procedures enable us to guarantee the reliability and consistency of our products. If you have further questions, please don't hesitate to contact us.



NanoIntegris Technologies, Inc.

c/o Raymor Industries

3765 La Verendrye

Boisbriand, Québec, Canada J7H 1R8

www.nanointegris.com

Telephone: 1-450-434-6266

Toll free: 1-866-650-0482

Toll free fax: 1-866-650-0482

Sales@nanointegris.com

NanoIntegris, a Raymor Company

T: 1-866-650-0482 | F: 1-514-434-9800

www.nanointegris.com | www.raymor.com | sales@nanointegris.com

2018 Product Catalog

Distributors



Sigma-Aldrich Co, LLC

- Website: www.sigmaaldrich.com
- Tel: 1-800-325-3010



Prime Machine Tools

- Website: www.primemachines.com
- Tel: +91-79-22901789
- Email: sales@primemachines.com



MKnano

- Website: www.mknano.com
- Tel: 416-509-4462
- Fax: 905-824-1259
- Email: sales@mknano.com



Chinwoo Tech. Corp.

- Website: www.chinwoo.co.kr
Contact: Hyounsub (Ian) Kim – Manager
- Tel: +82 31 777 1277 | Fax: +82 31 777 1288 | Mobile: +82 10 8976 5957
- Email: nanovea@chinwoo.co.kr



New Metals and Chemicals Corporation

- Website: www.newmetals.co.jp
- Tokyo - T. Itoh: ito@newmetals.co.jp
Tel: +81-3-5202-5624
- Osaka - K.
Kuriya: kuriya@newmetals.co.jp
Tel: +81-6-6202-5108

ANT Co.

- Website: www.ant-korea.com
Contact: Edward Cha.
- Tel: +82-2-555-6405
- Email: edward@ant-korea.com



WESi Technology Inc. (China and Taiwan)

- Website: www.wesitechnology.com
- Tel: +86-21- 5058-0051
Email: info@wesitechnology.com

KANTO CHEMICAL CO, INC

- Website: www.kanto.co.jp
- Email: reag-info@gms.kanto.co.jp
Tel: +81-03-3663-7631
- Contact: Shigeto Suga

XF Nano (China)

- Website: www.xfnano.com
- Tel: +86-25-6825-6969
- Email: sale@xfnano.com

OPTO SCIENCE, INC

- Website: www.optoscience.com
- Tokyo Office - M.
MURATA: murata@optoscience.com
Tel: +81-03-3356-1064
- Osaka Office - T.
OGAWA: ogawa@optoscience.com
Tel: +81-6-6305-2064

Boson Technologies (China)

- Website: www.BosonTech.com.cn
- Tel: +86-21-33520926
- Email: sales@bosontech.com.cn

NanoIntegris, a Raymor Company

T: 1-866-650-0482 | F: 1-514-434-9800

www.nanointegris.com | www.raymor.com | sales@nanointegris.com

RAYMOR
NanoIntegris
2018 Product Catalog

Table of Contents | Product List

Item SKU#	PRODUCT	Page
	IsoNanotubes-S: Semiconducting	
	Single-Wall Carbon Nanotubes (SWNT)	6
1101	99% Aqueous Solution	6
1102	99% Powder	6
1103	98% Aqueous Solution	6
1104	98% Powder	6
1105	95% Aqueous Solution	6
1106	95% Powder	6
1107	90% Aqueous Solution	6
1108	90% Powder	6
1117	(6,5) Chirally Enriched Solution	6
1130	(6,5) Chirally Enriched Powder	6
	IsoSol-S: Semiconducting	
	Single-Wall Carbon Nanotubes (SWNT)	8
1111	100% IsoSol-S100 Aromatic Solution	8
1112	100% IsoSol-S100 Powder	8
	IsoNanotubes-M: Metallic SWNTs	10
1201	99% Aqueous Solution	10
1202	99% Powder	10
1203	98% Aqueous Solution	10
1204	98% Powder	10
1205	95% Aqueous Solution	10
1206	95% Powder	10
1207	90% Aqueous Solution	10
1208	90% Powder	10
1209	70% Aqueous Solution	10
1210	70% Powder	10
	Unseparated Arc-Discharge SWNTs	12
1301	Super PureTubes Aqueous Solution	12
1302	Super PureTubes Powder	12
1303	PureTubes Aqueous Solution	12
1304	PureTubes Powder	12
1305	P2 Powder	12
1306	AP Powder	12

2018 Product Catalog

	Product Kits	14
1401	Starter Aqueous Solution	14
1402	Starter Powder	14
1403	Sample Aqueous Solution	14
1404	Sample Powder	14
1405	Premier Aqueous Solution	14
1406	Premier Powder	14
	HiPco Small Diameter SWCNTs	15
1601	Raw Fluffy Powder	15
1602	Raw Wet Cake	15
1603	Purified	15
1604	Super Purified	15
1613	HiPco(β) SWCNT - Raw	17
1614	HiPco(β) SWCNT – Purified	17
	HiPco Electronically Sorted Nanotubes	18
1605	Metallic Enriched Aqueous Solution	18
1606	Metallic Enriched Powder	18
1607	Semiconducting Enriched Aqueous Solution	18
1608	Semiconducting Enriched Powder	18
	PureWave Graphene: Crumpled Graphene Nanoplatelets	20
1708	PureWave Graphene Powder	
	PureSheets: 1-4 Layer Exfoliated Graphene Nanoplatelets	22
1701	MONO Aqueous Solution	22
1702	MONO Powder	22
1703	QUATTRO Aqueous Solution	22
1704	QUATTRO Powder	22
	PlasmaTubes SWNTs	24
1801	Super Purified (SPT-220) Aqueous Solution	24
1802	Super Purified (SPT-220) Powder	24
1803	Semi-Purified (RN-220) Aqueous Solution	24
1804	Semi-Purified (RN-220) Powder	24
1805	Air-Oxidized (RN-120) Aqueous Solution	24
1806	Air-Oxidized (RN-120) Powder	24
1807	Raw (RN-020) Aqueous Solution	24
1808	Raw (RN-020) Powder	24
1809	Raw (RN-000) Aqueous Solution	24
1810	Raw (RN-000) Powder	24

NanoIntegris, a Raymor Company

T: 1-866-650-0482 | F: 1-514-434-9800

www.nanointegris.com | www.raymor.com | sales@nanointegris.com

RAYMOR NanoIntegris

2018 Product Catalog

	Boron Nitride Nanotubes (BNNT)	26
1907	Purified BNNT	26
1905	Hexagonal Boron Nitride Flakes	26
1910	Multi-wall BNNT	27
	Transparent Conductive Ink	28
2004	Silver Nanowire Ink	28
	Double and Multi-Walled Carbon Nanotubes	29
2117	99%, Powder	29
2118	99%, Solution	29
2103	95%, <8nm OD, Powder	29
2105	95%, 10-20nm OD, Powder	29
2115	95%, 10-20nm OD, Solution	29
	Services	31
2201	Powder/ Thick Film Creation	31
2202	Dispersion: Aqueous	31
2203	Dispersion: Non-Aqueous Solution	31
2204	Small-scale Testing	31
2205	Production Scale-up	31
2206	Product Development/ R&D	31
2207	Scientific Consultation	31
2209	Solution Processing Fee	31
	Inks for Printed Electronics	32
1121	Xdi-dcs Dielectric Ink	32
1122	Xdi-dcs Ink & IsoSol-S100	32



RAYMOR
NanoIntegris
2018 Product Catalog


IsoNanotubes-S: Semiconducting SWCNT

Description: Semiconducting single-wall carbon nanotubes are a NanoIntegris exclusivity. They are sorted using a scalable production process. Offered in an aqueous solution or a powder.

Semiconducting carbon nanotubes may replace or complement traditional semiconductors in both high-performance and low-cost thin film transistor (TFT) devices. Today, TFTs are most commonly used in the backplanes of LCD and OLED displays. As the flexible electronics industry matures, TFTs will likely be incorporated into a much wider range of commercial electronics.

Please note that SKU# 1117 and 1118 use SWCNT from the CoMoCAT process.

Specifications:



RAYMOR NanoIntegris

2018 Product Catalog


Parameter	IsoNanotubes-S
Available Purities	90%, 95%, 98%, 99% and 99.9%
Metal content from TGA	< 1%
Length from AFM	0.3-4 μm
Solvent	Water or Toluene
Diameter range	1-1.7 nm

Prices:

SKU#	Product	Price (USD)	Quantity (mg)
1101	IsoNanotubes-S 99% Aqueous	800	1.0
1102	IsoNanotubes-S 99% Powder	850	1.0
1103	IsoNanotubes-S 98% Aqueous	450	1.0
1104	IsoNanotubes-S 98% Powder	500	1.0
1105	IsoNanotubes-S 95% Aqueous	300	1.0
1106	IsoNanotubes-S 95% Powder	350	1.0
1107	IsoNanotubes-S 90% Aqueous	200	1.0
1108	IsoNanotubes-S 90% Powder	250	1.0
1115	IsoNanotubes-S 99.9% Aqueous	1200	1.0
1116	IsoNanotubes-S 99.9% Powder	1250	1.0
1117	(6,5) Chirally Enriched 95% Aqueous	1150	200.0
1118	(6,5) Chirally Enriched 95% Powder	1150	1.0g

Significant discounts are available for large size orders. Please contact our sales team for special quote.

Did you know?





Semiconducting carbon nanotubes may replace or complement traditional semiconductors in both high-performance and low-cost thin film transistor (TFT) devices. Today, TFTs are most commonly used in the backplanes of LCD and OLED displays. As the flexible electronics industry matures, TFTs will likely be incorporated into a much wider range of commercial electronics.

IsoSol-S: Semiconducting SWCNT

Description: IsoSol-S100 is the award-winning semiconducting carbon nanotube ink with the highest purity ever measured coupled to a fully scalable manufacturing route that does not rely on ultracentrifugation. When processed optimally, our solution lead to carbon nanotube thin-film transistors (SWCNT-TFT) that deliver average mobilities of 10-20 cm²/V/s and current ON/OFF ratios of 10³-10⁶ (on Si/SiO₂) as well as average current densities of 1-10 μA/μm, sufficient in principle to drive organic light-emitting diodes (OLEDs). The semiconducting single-wall carbon nanotubes are sourced from the highly-scalable RF-plasma process, and separated using conjugated polymer extraction.

Specifications:



Property	Value
Optical Purity	>99.9%
Irkis Ratio (I)	>0.5
Phi Value (Φ)	>0.39
Nanotube Concentration	> 0.01mg/mL
Surfactant : Nanotube Concentration	<4
Standard Solvent Media	Toluene
Shelf Life	3 months

RAYMOR NanoIntegris


2018 Product Catalog

Prices:

SKU#	Product	Price (USD)	Quantity (mg)
1111	100% IsoSol-S100 Polymer Wrapped Powder	795	1.0
1112	100% IsoSol-S100 Aromatic Solution	695	1.0

Significant discounts are available for large size orders. Please contact our sales team for special quote.

Did you know?




Transistors can be printed using roll-to-roll (R2R) gravure on flexible polyethylene terephthalate (PET) substrate with silver gate electrodes and a high-k BaTiO₃ dielectric layer. The IsoSol-S100 was inkjet-printed at 50 mg/L along with silver source/ drain contacts, defining a channel of 1000 μm x 150 μm . The on/off ratios were determined to be 7×10^4 with a subthreshold swing of 1.8V/decade, a threshold voltage of -2.5V, and a mobility of 6cm²/Vs. **Reference:** ACS Appl. Mater. Interfaces 2016, 8, 27900-27910.

RAYMOR
NanoIntegris
2018 Product Catalog

IsoNanotubes-M: Metallic SWCNT

Description: Metallic single-wall carbon nanotubes are a NanoIntegris exclusivity. They are sorted using our patented density gradient ultracentrifugation (DGU) process and offered in an aqueous solution or a powder.

Specifications:



RAYMOR NanoIntegris

2018 Product Catalog


Parameter	IsoNanotubes-M
Available Purities	70%, 90%, 95%, 98% and 99%
Metal content from TGA	< 1%
Length range from AFM	0.3-4 μ m
Solvent	Water
Diameter range	1-1.7 nm

Prices:

SKU#	Product	Price (USD)	Quantity (mg)
1201	IsoNanotubes-M 99% Aqueous	900	1.0
1202	IsoNanotubes-M 99% Powder	950	1.0
1203	IsoNanotubes-M 98% Aqueous	700	1.0
1204	IsoNanotubes-M 98% Powder	750	1.0
1205	IsoNanotubes-M 95% Aqueous	400	1.0
1206	IsoNanotubes-M 95% Powder	450	1.0
1207	IsoNanotubes-M 90% Aqueous	300	1.0
1208	IsoNanotubes-M 90% Powder	350	1.0
1209	IsoNanotubes-M 70% Aqueous	200	1.0
1210	IsoNanotubes-M 70% Powder	250	1.0

Significant discounts are available for large size orders. Please contact our sales team for special quote.

Did you know?




Owing to their processability, stability, and high conductivity, carbon nanotubes have received significant attention from electronics-industry researchers over the past several years as an alternative to ITO. However, development work with nanotubes has been largely precluded by the unavoidable electronic polydispersity of as-grown CNTs. NanoIntegris has effectively solved this polydispersity problem—by separating as-grown nanotubes via DGU, we can produce large quantities of uniform metallic CNTs with up to 99% purity

Unseparated Arc-Discharge SWCNT

Description: These single-wall carbon nanotubes are purified to values of 95-98% (nanotube content measured via the Itkis ratio). PureTubes offer metal content lower than 5% and SuperPureTubes offer less than 2% metal content (as per TGA). These are our highest purity un-sorted nanotubes. Sourced from the arc-discharge processes and offered in solution or dry powder.

Specifications:



NanoIntegris, a Raymor Company

T: 1-866-650-0482 | F: 1-514-434-9800

www.nanointegris.com | www.raymor.com | sales@nanointegris.com

RAYMOR NanoIntegris

2018 Product Catalog


Parameter	PureTubes	SuperPureTubes
G/D ratio from Raman	> 30	> 30
Metal content from TGA	< 4%	1%
Length range from AFM	0.1-4 μm	0.1-4 μm
Carbon impurities	<5%	<5%
Diameter range	1-1.7 nm	1-1.7 nm
Itkis index from Optical Absorbance	> 0.2	> 0.2

Prices:

SKU#	Product	Price (USD)	Quantity (mg)
1301	Super PureTubes Aqueous	2500	1.0
1302	Super PureTubes Powder	2600	1.0
1303	PureTubes Aqueous	960	1.0
1304	PureTubes Powder	1060	1.0
1305	P2 Powder	350	1.0g
1306	AP Powder	135	3.0g
1307	COOH-Functionalized P2	1000	1.0g
1308	COOH-Functionalized AP	750	1.0g

Significant discounts are available for large size orders. Please contact our sales team for special quote.

Did you know?



Photonic devices are widely used for optical communications, spectroscopy, and precision surgery (e.g. medical lasing). Materials which exhibit strong nonlinear electro-optical behaviors are required for most photonic applications. Ideally, these materials should exhibit fast response times, absorb over a broad wavelength range, and exhibit low optical loss. Nanotubes are one of a handful of materials in existence which satisfy these property requirements.

2018 Product Catalog

Product Kits

Description: Dozens of researchers [publish groundbreaking research using our tubes](#) every year (over 500 publications as of January 2016). Make sure the next major paper published in Nature, JACS, ACS Nano or Nano Letters is yours! If you are starting a new project and you are not sure which nanotube product will perform best, try a starter kit, sample or premium kit and make sure to save over 60% of the retail price if items were sold separately.

Starter kits include: **2 mg of IsoNanotubes 90%-M**
 2 mg of IsoNanotubes 90% - S
 50 mg of PureTubes

Sample kits include: **2 mg of IsoNanotubes 70%-M, 2 mg of IsoNanotubes 90%-M**
 2 mg of IsoNanotubes 90% - S, 2 mg of IsoNanotubes 95%-S
 100 mg of PureTubes

Premier kits include **1 mg of IsoNanotubes 99%-M, 1 mg of IsoNanotubes 98%-M**
 1 mg of IsoNanotubes 99% - S, 1 mg of IsoNanotubes 98%-S
 50 mg of PureTubes.

Specifications: see section above for specifications of IsoNanotubes – M, IsoNanotubes – S and PureTubes.



Prices:

SKU#	Product	Price (USD)	Quantity (mg)
1401	Starter Aqueous Solution	900	See above
1402	Starter Powder	950	See above
1403	Sample Aqueous Solution	700	See above
1404	Sample Powder	750	See above
1405	Premier Aqueous Solution	400	See above
1406	Premier Powder	450	See above

NanoIntegris, a Raymor Company


T: 1-866-650-0482 | F: 1-514-434-9800

www.nanointegris.com | www.raymor.com | sales@nanointegris.com

HiPco Small Diameter SWCNT

Description: HiPco SWCNT represent a benchmark for small diameter nanotubes both in the academic community and for industrial and commercial applications. They are synthesized in the original fluidized bed process developed by Nobel laureate Professor Richard Smalley and located at Rice University. We offer three grades of purity for this unique material. Raw HiPco SWCNT contain up to 15% of iron, whereas Purified HiPco SWCNT and Superpurified HiPco SWCNT contain less than 10 and 5% iron respectively.

Specifications:



NanoIntegris, a Raymor Company

T: 1-866-650-0482 | F: 1-514-434-9800

www.nanointegris.com | www.raymor.com | sales@nanointegris.com

RAYMOR NanoIntegris

2018 Product Catalog


Parameter	Raw	Purified	SuperPurified
Residual Iron Content	<35%	<15%	<5%
Length	0.1-1 μm	0.1-1 μm	0.1-1 μm
Diameter	0.8-1.2 nm	0.8-1.2 nm	0.8-1.2 nm

Prices:

SKU#	Product	Price (USD)	Quantity (g)
1601	Raw Fluffy Powder	550	1.0
1602	Raw Wet Cake	550	1.0
1603	Purified	800	1.0
1604	Super Purified	2200	1.0

Discounts are available for large size orders. Please contact our sales team for special quote.

Did you know?




Chemical sensors are used for many purposes, such as environmental hazard screening, explosives detection, product characterization, and medical testing. The electronic properties of SWNTs can change significantly when gases and bio-molecules are adsorbed to their surface. These changes can be detected in resistor, transistor, or capacitor devices. A principle advantage of TFT SWNT sensors in particular is that they respond to analyte surface coverage, as opposed to conventional sensors, which respond to analyte concentration.

HiPco(β): Small Diameter SWCNT

Description: NanoIntegris is very proud to offer a new grade of HiPco SWCNT. These small diameter single-wall carbon nanotubes are synthesized using the high-pressure carbon monoxide process developed at Rice University. A new reactor, recently launched, produces the same material, but at a different location and at larger scale.

Specifications:



Metric	Value	Method of Analysis
Diameter	0.6 – 1.1 nm	Raman, TEM
Length	400 – 700 nm	TEM
Color	Black	-
Ash content of Raw SWCNT	< 10 wt%	TGA
Ash content of Purified SWCNT	< 1 wt%	TGA
Carbon Purity	> 90%	TGA
Average G/D Ratio	~ 35	Raman (514 nm)

Prices:

SKU#	Product	Price (USD)	Quantity (g)
1613	HiPco(β) SWCNT - Raw	299	1.0
1614	HiPco(β) SWCNT - Purified	549	1.0

NanoIntegris, a Raymor Company


T: 1-866-650-0482 | F: 1-514-434-9800

www.nanointegris.com | www.raymor.com | sales@nanointegris.com

HiPco Semiconducting/Metallic Sorted SWCNT

Description: HiPco SWCNT represent a benchmark for small diameter nanotubes both in the academic community and for industrial and commercial applications. They are synthesized using the fluidized bed process developed by Nobel laureate Richard Smalley. We offer three grades of purity for this unique material. Raw HiPco SWCNT contain up to 15% of iron, whereas Purified HiPco SWCNT and Superpurified HiPco SWCNT contain less than 10 and 5% iron respectively.

Specifications:




Product	Semiconducting Enrichment	Metallic Nanotube Content	Diameter	Residual Fe Catalyst
HiPco Semiconducting SWCNT	>98%	<2%	0.8-1.2nm	<5% by weight
HiPco Metallic SWCNT	<5%	>95%	0.8-1.2nm	<5% by weight

NanoIntegris, a Raymor Company

T: 1-866-650-0482 | F: 1-514-434-9800

www.nanointegris.com | www.raymor.com | sales@nanointegris.com

RAYMOR
NanoIntegris
2018 Product Catalog



Prices:

SKU#	Product	Price (USD)	Quantity (g)
1605	Metallic Enriched Aqueous Solution	600	1.0
1606	Metallic Enriched Powder	650	1.0
1607	Semiconducting Enriched Aqueous Solution	500	1.0
1608	Semiconducting Enriched Powder	550	1.0

2018 Product Catalog

PureWave Graphene

Description: Our PureWave Graphene Nanoplatelets (GNP) material is composed of thin, highly dispersible graphene nanoplatelets with very low oxygen content. The turbostratic and wavy morphology of the material leads to an unequalled ability to be dispersed in a variety of solvents and resins. Furthermore, the unique plasma process used to grow our PureWave Graphene Nanoplatelets is easy to scale and thus produces a low cost product. Our unique and protected growth process based on plasma enables us to produce this material at over 200 g/hour, enabling commercial graphene applications such as energy storage and conductive inks as well as elastomers and resins.

Their non-planar morphology imparts these GNPs with a comparatively weak tendency to aggregate relative to planar exfoliated GNPs. This and their small size results in excellent dispersibility in solvents and polymer matrices even though they are pristine. This characteristic coupled with their excellent electrical conductivity makes PureWave graphene especially effective in conductive composites, coatings, inks and energy storage devices such as lithium ion battery anodes, cathodes and in supercapacitors. In several important polymeric matrices and even in solvent, PureWave graphene has been found to give conductive formulations (percolate) at lower concentrations than exfoliated GNPs. Addition of PureWave plasma GNPs to formulations containing larger planar exfoliated GNPs gives a substantial synergistic boost to conductivity, indicating that the small PureWave GNPs are very effective at electrically connecting the larger particles. Likewise, plasma GNPs are especially effective as the conductive component in cathode and anode coatings.

Further reading: Types and Applications of Graphene: Introducing PureWave Plasma GNPs. N. Vanier et al. (2014) <https://www.printedelectronicsworld.com/articles/10940/types-and-applications-of-graphene-introducing-purewave-plasma-gnps>

Specifications:

Physical Characterization	
Parameter	Value
G/D ratio from Raman (514 nm)	3
2D/G ratio from Raman	0.8-1
BET Specific Surface Area	> 400 m ² /g
Average Number of Layers	6
Flake Size	150-200 nm
Oxygen Content	1%
Average Flake Thickness	2.4 nm
Carbon Content (TGA)	> 98%
Metal Impurities (ICP-MS)	< 300 ppm
Performance in Formulation (screen printed)	
Parameter	Value
Resistivity (Ω/□/mil)	1.0 - 5.0


NanoIntegris, a Raymor Company

T: 1-866-650-0482 | F: 1-514-434-9800

www.nanointegris.com | www.raymor.com | sales@nanointegris.com

RAYMOR NanoIntegris

2018 Product Catalog



A: TGA and TGA time derivative of PureWave Graphene Nanoplatelets.
B: Raman spectra acquired at 514 nm excitation wavelength. **C:** TEM microscopy image of Graphene Nanoplatelets.

Prices:

SKU#	Product	Price (USD)	Quantity (g)
1708	PureWave Graphene Nanoplatelets (GNP)	100	2.0
1708	PureWave Graphene Nanoplatelets (GNP)	500	20.0
1708	PureWave Graphene Nanoplatelets (GNP)	1,500	100.0

RAYMOR NanoIntegris

2018 Product Catalog


PureSheets: 1-4 Layer Graphene Nanoplatelets

Description: PureSheets graphene nanoplatelets (GNP) were mechanically exfoliated from graphite. Our graphene solutions are post-processed to remove thicker platelets that were not properly exfoliated. This step ensures the optimal dispersion and stability of our products. PureSheets-MONO GNP contain, on average, 1-3 graphene layers according to AFM while PureSheets-QUATTRO contain, on average, 4-6 layers of graphene.

Specifications:



- MONO Grade
- QUATTRO Grade



NanoIntegris, a Raymor Company

T: 1-866-650-0482 | F: 1-514-434-9800

www.nanointegris.com | www.raymor.com | sales@nanointegris.com

RAYMOR NanoIntegris

2018 Product Catalog


Parameter	MONO	QUATTRO	Measurement
Single Layer Content	27%	6%	AFM
Double Layer Content	48%	23%	AFM
Triple Layer Content	20%	27%	AFM
4+ Layer Content	5%	44%	AFM
Average Flake Area	10000 nm ²	10000 nm ²	AFM
Solution Type	Acqueous	Acqueous	n/a
Graphene Concentration	0.05 mg/ml	0.05 mg/ml	n/a
Surfactant Concentration	2% w/v	2% w/v	n/a

Prices:

SKU#	Product	Price (USD)	Quantity (mg)
1701	MONO Aqueous	625	400.0
1702	MONO Powder	675	400.0
1703	QUATTRO Aqueous	400	400.0
1704	QUATTRO Powder	450	400.0
1708	PureWave graphene nanoplatelets	100	5.0g

Significant discounts are available for large size orders. Please contact our sales team for special quote.

Did you know?



Large-diameter semiconducting SWNTs are good absorbers and emitters of light in the infrared. Moreover, high-purity SWNT thin films have been demonstrated to be photoconductive and photo-luminescent under NIR illumination.


IR sensors/emitters are useful for a number of military and civilian applications.

2018 Product Catalog

PlasmaTubes SWCNT

Description: Using a patented plasma torch process, Raymor Nanotech produces raw single-wall carbon nanotubes (SWCNT) at high rates, enabling the lowest prices on the market. As shown by the detailed analysis below, plasma-grown SWCNT display a high graphitization level, diameters (0.9-1.5 nm) and lengths (0.3-4 μm) close to those of laser- and arc-grown SWCNT. The purity of the raw SWCNT (RN-020) is comparable to the purity of the best arc-discharge SWCNT on the market. The following sections will display the information gathered by thermogravimetric analysis, Raman spectroscopy and optical absorption on the raw (RN-020), purified (RN-220) and SuperPurified (SPT-220) SWCNT. The Super Purified SWCNT (SPT-220) product has a nanotube purity of 95-99%. This material is provided in an aqueous surfactant solution with a nanotube concentration of 0.50 mg/ml or a surfactant-removed thick film. Regarding the difference between RN-020, and RN-000 raw SWCNT, the first are grown using a Co-Fe-Ni catalyst mixture whereas the latter are grown using a Ni-Y catalyst mixture. Unless requested, RN-020 is used for all purified and separated grades.

Specifications:



Parameter	Measurement for RN-020 & RN-000	Typical range
G/D ratio with BWF subtraction	Raman spectroscopy at 514 nm	> 35
G/D ratio without BWF subtraction	Raman spectroscopy at 514 nm	> 30
Ash content	Thermogravimetric analysis	27%
1st oxidation peak	Thermogravimetric analysis	400 °C
2nd oxidation peak	Thermogravimetric analysis	690 °C
I _t k _i s index	Optical absorption	0.06-0.08
Carbon Purity	Thermogravimetric analysis	80-85%

NanoIntegris, a Raymor Company

T: 1-866-650-0482 | F: 1-514-434-9800

www.nanointegris.com | www.raymor.com | sales@nanointegris.com

RAYMOR
NanoIntegris
2018 Product Catalog

Parameter	Measurement for RN-220	Typical range
G/D ratio with BWF subtraction	Raman spectroscopy at 514 nm	> 80
G/D ratio without BWF subtraction	Raman spectroscopy at 514 nm	> 70
Ash content	Thermogravimetric analysis	21%
1st oxidation peak	Thermogravimetric analysis	580 °C
2 nd oxidation peak	Thermogravimetric analysis	650 °C
Itkis index	Optical absorption	0.08-0.1

Parameter	Measurement for SPT-220	Typical range
G/D ratio with BWF subtraction	Raman spectroscopy at 514 nm	> 40
G/D ratio without BWF subtraction	Raman spectroscopy at 514 nm	> 35
Ash content	Thermogravimetric analysis	1-3%
1st oxidation peak	Thermogravimetric analysis	580 °C
2 nd oxidation peak	Thermogravimetric analysis	650 °C
Itkis index	Optical absorption	> 0.2

Prices:

SKU#	Product	Price (USD)	Quantity (g)
1801	Super Purified (SPT-220) Aqueous	1000	1.0
1802	Super Purified (SPT-220) Powder	1100	1.0
1803	Semi-Purified (RN-220) Aqueous	155	1.0
1804	Semi-Purified (RN-220) Powder	85	1.0
1805	Air-Oxidized (RN-120) Aqueous	100	1.0
1806	Air-Oxidized (RN-120) Powder	42	1.0
1807	Raw (RN-020) Aqueous	35	1.0
1808	Raw (RN-020) Powder	20	1.0
1809	Raw (RN-000) Aqueous	35	1.0
1810	Raw (RN-000) Powder	20	1.0



Significant discounts are available for large size orders. Please contact our sales team for special quote.

RAYMOR
NanoIntegris
2018 Product Catalog

Purified Few-Walled Boron Nitride Nanotubes (BNNT)

Description: Boron Nitride Nanotubes (BNNT) are a brand new addition to our catalog. Their extreme resistance to high temperatures (up to 900°C in air) and optical transparency make them an ideal filler for the structural reinforcement of glass, ceramics and metal matrices. Furthermore, their high thermal conductivity may enable BNNT to solve one of the main challenges in electronic engineering: thermal management of ever smaller circuit components. Our BNNT are synthesized using a patented plasma torch process, leading to scalable production.

Specifications:



Attribute	Measurement
BNNT Purity	> 69 % by weight
Number of Walls and Diameter	1 - 5 walls typical / ~ 5 nm diameter
Tube Length	≤ 200 μm
Surface Area	100-1800 m ² /g
Bundles	≤ 5 tubes across
Electrical Conductivity / Band Gap	Large band gap insulator ~5.5 eV
Total h-BN & BNH content	≤ 20 %
Tap Density	~0.25mg/cm ³
Thermal Stability	Stable in highly oxidizing environments, e.g. ~900 °C in air.

Prices:

SKU#	Product	Price (USD/g)	Quantity (g)
1907	Purified BNNT Powder	1650	0.5-49
1907	Purified BNNT Powder	1350	49-99
1907	Purified BNNT Powder	1065	100+
1905	Hexagonal Boron Nitride Flakes	1	100.0

NanoIntegris, a Raymor Company
T: 1-866-650-0482 | F: 1-514-434-9800
www.nanointegris.com | www.raymor.com | sales@nanointegris.com


RAYMOR NanoIntegris

2018 Product Catalog

Multi-Walled Boron Nitride Nanotubes (MW-BNNT)

Description: Multi-Walled Boron Nitride Nanotubes (MW-BNNT) may provide the answer to your microelectronic circuit overheating problems! Their extreme resistance to high temperatures (up to 900°C in air) and optical transparency make them ideal for the structural reinforcement of glass, ceramics and metal matrices. In a thermal conductivity epoxy composite comprising 0.5wt% MW-BNNT, thermal conductivity was enhanced twofold (from 1.3 to 3.4 W/mK) compared to an epoxy contained h-BN only. It is projected that this material could provide a thermal conductivity up to 10 W/mK in such epoxy composites. Finally, the electric breakdown voltage of these composites was measured at 26.3 kV/mm. Our MW-BNNT are synthesized by a CVD process.

Specifications:



Metric	Value
BET Surface Area	> 30 m ² /g
Diameter	50-100 nm
Length	10 - 20 micron

Prices:

SKU#	Product	Price (USD/g)	Quantity (g)
1910	BNNT-MW Powder	625	1 - 99
1910	BNNT-MW Powder	500	100 - 199
1910	BNNT-MW Powder	300	200 -299

RAYMOR NanoIntegris

2018 Product Catalog

Silver Nanowire Ink for Transparent Conductive Film

Description: We offer a new high performance silver nanowire ink for transparent conductive film applications. Our highly conductive transparent ink works can be coated with slot die coating and micro-gravure coating, along with other application techniques.

Specifications:



Performance Index	Testing method	Testing Result
AgNW diameter / length	TEM / SEM	30nm / 20µm
Appearance	Visual	Gray suspension
AgNW content (wt%)	Wet Combustion Method	0.2 - 0.3
Density (g/ml)	Densimeter	1.05
Viscosity (cps) @ 25°C	Rotor rotational viscometer	5 -30 cps
Curing temperature	130°C for 3-5 min	
Sheet resistance ($\Omega/\square$)	Four Probe Method	50 - 100
Transmittance (%)	WGW	90 - 91
Haze (%)	WGW	1 - 1.5
Surface hardness (H)	Pencil hardness testing device	2 - 3H
Adhesion	3M 600 tape, pull vertically	Not shed

Prices:

SKU#	Product	Price (USD)	Quantity
2004	Silver Nanowire Ink	350	1 g / 100 ml

Significant discounts are available for large size orders. Please contact our sales team for special quote.

NanoIntegris, a Raymor Company
T: 1-866-650-0482 | F: 1-514-434-9800
www.nanointegris.com | www.raymor.com | sales@nanointegris.com

2018 Product Catalog

Double and Multi-Walled Carbon Nanotubes (MWCNT)

Description: Multi-walled carbon nanotubes (MWCNT) are the ideal filler for the structural reinforcement or improvement in electrical conductivity of polymers, elastomers and epoxy. These CVD-grown nanotubes display high purity and low metal content, a controlled diameter range, good length and aspect ratio but most importantly are priced to suit both research and industrial needs.

Specifications:



Parameter	99% Purity MWNT	<8 nm MWNT	10-20 nm MWNT
MWNTs Outer Diameter	<20nm	<8 nm	10-20 nm
MWNTs Inside Diameter	4nm	2-5 nm	3-5 nm
MWNTs Ash	0 wt%	<1.5 wt%	<1.5 wt%
MWNTs Purity	>99 wt%	>95 wt%	>95 wt%
WNTs Length	1-12 μ m	10-30 μ m	10-30 μ m
Source Material	CVD	CVD	CVD
MWNTs Specific Surface Area	n/a	500 m ² /g	233 m ² /g
MWNTs Electrical Conductivity	n/a	>100 S/cm	>100 S/cm
MWNTs Bulk density	n/a	0.27 g/cm ³	0.22 g/cm ³
MWNTs True density	n/a	~2.1 g/cm ³	~2.1 g/cm ³

NanoIntegris, a Raymor Company

T: 1-866-650-0482 | F: 1-514-434-9800

www.nanointegris.com | www.raymor.com | sales@nanointegris.com

RAYMOR NanoIntegris


2018 Product Catalog

Prices:

SKU#	Product	Price (USD)	Quantity (g)
2117	99%, Powder	125	5
2118	99%, Solution	175	1
2103	95%, <8nm OD Powder	65	5
2103	95%, <8nm OD Solution	360	100
2105	95%, 10-20nm OD Powder	50	5
2115	95%, 10-20nm OD Solution	100	1

Significant discounts are available for large size orders. Please contact our sales team for special quote.

Did you know?



CNTs have proven useful for targeted drug delivery. Anti-cancer drugs may be delivered more efficaciously and with fewer systemic side-effects using a "smart" nanotechnology platform than by conventional methods. Carbon nanotubes represent one such promising platform, due to their strong absorbance in the so-called therapeutic infrared window (between 700-1100 nm, depending on body tissue type).

2018 Product Catalog

Services

Description: Take advantage of our 10 years of expertise in the nanotechnology and specifically in the nanotubes industry. Raymor and NanoIntegris constitute one of the world's largest and most trusted provider of single-walled carbon nanotubes. We have produced, purified, separated, functionalized and dispersed nanotubes using the best methods in a variety of media and solvents. Please contact us for any special projects or inquiries. It will be our pleasure to provide you with the simplest and most cost effective solution. The following list is subset of available services so please do not hesitate to contact us for your special requirements.



Prices:

SKU#	Product	Price (USD)
2201	Powder/ Thick Film Creation	50
2202	Dispersion: Aqueous	Ask for quote
2203	Dispersion: Non-Aqueous Solution	Ask for quote
2204	Small-scale Testing	Ask for quote
2205	Production Scale-up	Ask for quote
2206	Product Development/ R&D	Ask for quote
2207	Scientific Consultation	Ask for quote
2209	Solution Processing Fee	Ask for quote

NanoIntegris, a Raymor Company

T: 1-866-650-0482 | F: 1-514-434-9800

www.nanointegris.com | www.raymor.com | sales@nanointegris.com

RAYMOR
NanoIntegris
2018 Product Catalog


Inks for Printed Electronics

Description: NanoIntegris has begun selling a complementary thin film transistor package. The materials package combines a new **dielectric ink (xdi-dcs)** developed at the XRCC with NanoIntegris' high purity, single-walled carbon nanotube ink (IsoSol-S100), to improve the overall performance of printed high-mobility p-type transistors. Compatibility between semiconducting and dielectric materials is critical for reliable processing and device performance.

One of the challenges that have limited the implementation of single-walled carbon nanotube based thin film transistors is that they exhibit considerable hysteresis with a non-zero threshold voltage when exposed/operated in air ambient. This new materials package addresses performance issues when the **Xdi-dcs ink** is used as a **dielectric**, and when further used as an encapsulant. Key to addressing issues such as hysteresis is the hydrophobic surface of the dielectric layer, which can eliminate water and mitigate charge trapping. Both materials have low viscosity, enabling spin coating, inkjet printing or aerosol deposition, and can be cured at or below 150 °C, making them compatible with the most popular low temperature substrate polyethylene terephthalate (PET).

Further details: J. Lefebvre et al. Hysteresis free carbon nanotube thin film transistors comprising hydrophobic dielectrics. Appl. Phys. Lett. 107, 243301 (2015).

Specifications:



RAYMOR
NanoIntegris
2018 Product Catalog

Solution Characteristics		
Metric	IsoSol Semiconductor	Xdi-dcs Dielectric
Viscosity at 25 °C	18-22 cps	8-9 cps
Surface Tension	25-30 mN/m	24-25 mN/m
Solvent System	organic	organic
Cure (thermal, air)	150 °C / 10 min	140 °C / 30 min
Shelf Life	up to 3 months	up to 2 months
Preferred Deposition Method	Aerosol	Spin-coating
Thin Film Transistor Characteristics		
Metric	Value (in air)	Value (encapsulated)
Mobility	40 cm ² /V/s	> 6 cm ² /V/s
On/Off Ratio	1 x 10 ⁴	1 x 10 ²
Threshold Voltage	4 ± 1 V	0 ± 1 V
Hysteresis	0.004 ± 0.03 V	0.004 ± 0.03 V
Dielectric Constant	4	4
Dielectric Layer Thickness	~ 500 nm	530 nm

Prices:

SKU#	Product	Price (USD)	Quantity (g)
1121	Xdi-dcs Dielectric Ink	16 \$/g	< 2000
1121	Xdi-dcs Dielectric Ink	12 \$/g	> 2000
1122	10g XDI-DCS Ink & 1mg IsoSol-S100	775	10 g Xdi-dcs+ 1 mg IsoSol

Please view page 8 for the details of IsoSol-S100 semiconducting ink.